GILUDUR





Synthetic Hard Plaster

EN ISO 6873 -3

For master- and working-models and as flask material in the plastic prosthetic technique, also for occluding dentition.

GILUDUR offers

- models of high strengthbalanced expansion
- consistant quality

Physical data

Setting time app. 10 min.
Compressive strength
after 24 hours min. 30 N/mm²
Hardness
after 24 hours min. 80 N/mm²
Setting expansion
after 2 hours max. 0.2 %

Working method

Mixing ratio
GILUDUR: water
Water: powder
Working temperature

Mixing ratio
100 g: 30 ml
0.3 ml/g
app. 23° C

Strew the powder in the given ratio into the water within 10 sec. and allow the mix to soak for 20 sec.

Mechanical mixing	30 sec
Hand mixing	60 sec

The initial consistency becomes more fluid during mixing and reaches optimum fluidity, when pouring the moulds on the vibrator.

Pouring time app. 6 min

Do not vibrate after this time so as not to disturb the crystallization.

Availability

Blue or yellow; special colours can be delivered, minimum batch quantity per shipment is 2 tons.

Packaging

25 kg
45 kg
20 kg
4,5 kg
18 kg

Shelf life

At least 1 year in well closed, vapour-proof packages.

Important

GILUDUR - like all plaster products - is sensitive to humidity and contact with the air should be avoided. If repackaged use only vapour-proof material e.g.

- plastic aluminium foil
- low-pressure-polythene containers with a wall thickness of min. 0.5 mm.

Container should be well sealed between usage and should be stored in a dry area.

GILUDUR which has been stored at extreme temperatures should be equilibrated at R.T. for several hours before use.

It is recommended that a model which has dried out during storage be soaked in water for approx. 2 min. before sawing or removing wax (injection method).

The recommendations are given to the best of our knowledge after careful control. We grant the quality of our products according to our specification. Any further liability cannot be accepted since the proper application of our products is outside of our control.

